



### **Product description**

Order number	65.163
Product description	Screw cap, material: HD-PE, natural, for tubes 75 x 13 mm, 1,000 piece(s)/bag

### **Product characteristics**

Suitable for	tubes 75 x 13 mm
--------------	------------------

### **Material & colours**

Product material	High Density Polyethylene (HD-PE)
Colour of product	natural

### **Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Packaging**

Minimum order qty.	5000
Type of smallest subpackaging	bag
Piece(s) / inner box	1000
Piece (s) / outer case	5000
Piece(s) / pallet	210000
Depth of case	398 mm
Width of case	398 mm
Height of case	250 mm
Case volume	0.0396 cbm
Weight of product	0.0014 kg
Weight of case	7.45 kg
EAN of inner box	4038917272223
EAN case	4038917272216

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

