



Product Specification
Safety-Multifly® BKF set 23G, 23G x 3/4", blue,
tube length: 200 mm, 1 piece(s)/blister



Product description

Order number	51.1640.945
Product description	Safety-Multifly® BKF set 23G, needle (ØxL): 23G x 3/4", coding colour: blue, tube length: 200 mm, with blood culture adapter set, set components: 1 Safety-Multifly® needle 23G with long tubing, DEHP-free (REF 85.1640.235) and 1 Blood Culture Adapter Universal (REF 14.1209), 1 piece(s)/blister, sterile, 1 piece(s)/case

Product characteristics

Design	blood culture adapter set
Coding colour	blue
Safety product	yes
Connection	blood culture bottle

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Safety-Multifly® BKF set 23G, 23G x 3/4", blue,
tube length: 200 mm, 1 piece(s)/blister

Page 2

Size

Diameter of needle (external)	23 G / 0,6 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	200 mm

Purity & certification

Product category	Medical device
Purity standard	sterile
Batched	yes

Packaging

Minimum order qty.	50
Type of smallest subpackaging	blister
Piece(s) / inner box	1
Piece (s) / outer case	50
Piece(s) / pallet	2100
Depth of case	388 mm
Width of case	378 mm
Height of case	235 mm
Case volume	0.0345 cbm
Weight of product	0.0118 kg
Weight of case	1.42 kg
EAN of inner box	4038917357500
EAN case	4038917357494

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-06-30

This document was prepared by EDP support and is valid without signature.



SARSTEDT