



Product description

Order number	85.1637.203
Product description	Safety-Multifly® needle, needle (ØxL): 20G x 3/4", coding colour: yellow, tube length: 80 mm, with Luer (f), DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case
PZN number	6797488

Product characteristics

Design	Luer (f)
Coding colour	yellow
Safety product	yes
Connection	Luer (m)

Size

Diameter of needle (external)	20 G / 0,9 mm
Diameter of tube	2.3 mm
Length of needle	19 mm
Length of tube	80 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Safety-Multifly® needle, 20G x 3/4", yellow, tube
length: 80 mm, 1 piece(s)/blister

Purity & certification

Product category	Medical device
CE certified	CE 0124 (DEKRA)
Purity standard	sterile, non-pyrogenic/endotoxin-free
Sterilisation	Gas sterilisation
Pyrogen-free	yes
Plasticiser	DEHP-free
Batched	yes

Packaging

Minimum order qty.	480
Type of smallest subpackaging	blister
Piece(s) / inner box	120
Piece (s) / outer case	480
Piece(s) / pallet	20640
Depth of box	224 mm
Width of box	214 mm
Height of box	171 mm
Depth of case	442 mm
Width of case	358 mm
Height of case	240 mm
Case volume	0.038 cbm
Weight of product	0.0023 kg
Weight of case	2.86 kg
EAN of inner box	4038917423854
EAN case	4038917423564

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

