



Product Specification
**Multifly® needle, 21G x 3/4", green, tube length:
80 mm, 1 piece(s)/blister**

Page 1

Product description

| | |
|---------------------|---|
| Order number | 85.1638.005 |
| Product description | Multifly® needle, needle (ØxL): 21G x 3/4", coding colour: green, tube length: 80 mm, with Multi-adapter, DEHP-free, 1 piece(s)/blister, sterile, non-pyrogenic/endotoxin-free, 120 piece(s)/case |

Product characteristics

| | |
|----------------|--------------------|
| Design | Multi-adapter |
| Coding colour | green |
| Safety product | no |
| Connection | membrane screw cap |

Size

| | |
|-------------------------------|---------------|
| Diameter of needle (external) | 21 G / 0,8 mm |
| Diameter of tube | 2.3 mm |
| Length of needle | 19 mm |
| Length of tube | 80 mm |

Purity & certification

| | |
|------------------|---------------------------------------|
| Product category | Medical device |
| CE certified | CE 0124 (DEKRA) |
| Purity standard | sterile, non-pyrogenic/endotoxin-free |
| Sterilisation | Gas sterilisation |
| Pyrogen-free | yes |
| Plasticiser | DEHP-free |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------|
| Minimum order qty. | 480 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 120 |
| Piece (s) / outer case | 480 |
| Piece(s) / pallet | 20640 |
| Depth of box | 224 mm |
| Width of box | 214 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-09-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
**Multifly® needle, 21G x 3/4", green, tube length:
80 mm, 1 piece(s)/blister**

| | |
|-------------------|---------------|
| Height of box | 171 mm |
| Depth of case | 442 mm |
| Width of case | 358 mm |
| Height of case | 240 mm |
| Case volume | 0.038 cbm |
| Weight of product | 0.0027 kg |
| Weight of case | 3.19 kg |
| EAN of inner box | 4038917423762 |
| EAN case | 4038917423472 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

