




**Product Specification**  
**Sample tube, Serum CAT, 10 ml, cap white, (LxØ):**  
**95 x 16.8 mm, with print**

**Page 1**



**Product description**

Order number	26.367
Product description	Sample tube, nominal volume: 10 ml, cap white, preparation: Clotting activator, push cap, (LxØ): 95 x 16.8 mm, round base, with print, 100 piece(s)/bag

**Product characteristics**

Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with print
Colour of print/label	white
Graduation	yes
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-01

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Sample tube, Serum CAT, 10 ml, cap white, (LxØ):**  
**95 x 16.8 mm, with print**

**Size**

Sample volume	10 ml
Diameter	16.8 mm
Length of product	95 mm
Length excluding cap	95 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	white
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0066 kg
Weight of case	7.27 kg
EAN of inner box	4038917081337
EAN case	4038917000475

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

