



## Product Specification

Sample tube, Citrate 9NC 0.106 mol/l 3.2%, 5 ml,  
cap green, (LxØ): 75 x 13 mm, with print

Page 1



### Product description

Order number	32.525.326
Product description	Sample tube, nominal volume: 5 ml, cap green, preparation: Citrate 9NC, push cap, (LxØ): 75 x 13 mm, round base, with print, 100 piece(s)/polystyrene container

### Product characteristics

Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	EU
Label/ Print	with print
Colour of print/label	green
Graduation	fill mark
Closure type	push cap
Base shape	round base
Recommended centrifugation	1800 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Sample tube, Citrate 9NC 0.106 mol/l 3.2%, 5 ml,  
cap green, (LxØ): 75 x 13 mm, with print**

**Page 2**

**Size**

Sample volume	5 ml
Diameter	13 mm
Length of product	75 mm
Length excluding cap	75 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	green
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	polystyrene container
Piece(s) / inner box	100
Piece (s) / outer case	500
Piece(s) / pallet	37500
Depth of case	458 mm
Width of case	348 mm
Height of case	115 mm
Case volume	0.0183 cbm
Weight of product	0.0038 kg
Weight of case	2.42 kg
EAN of inner box	4038917109277
EAN case	4038917012454

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**