




Product Specification
Sample tube, EDTA K3E, 4 ml, cap red, (LxØ): 75 x 12 mm, with print

Page 1



Product description

Order number	36.338
Product description	Sample tube, nominal volume: 4 ml, cap red, preparation: K3 EDTA, push cap, (LxØ): 75 x 12 mm, round base, with print, 100 piece(s)/bag

Product characteristics

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with print
Colour of print/label	red
Graduation	yes
Cap	Skirted push cap
Closure type	push cap
Base shape	round base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
Sample tube, EDTA K3E, 4 ml, cap red, (LxØ): 75 x 12 mm, with print

Page 2

Size

Sample volume	4 ml
Diameter	12 mm
Length of product	75 mm
Length excluding cap	75 mm

Material & colours

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	red
Tube material	Polypropylene (PP)

Purity & certification

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

Packaging

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0033 kg
Weight of case	7.32 kg
EAN of inner box	4038917081238
EAN case	4038917000376

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT