




**Product Specification**  
**Sample tube, Fluoride/heparin FH, 2 ml, cap yellow, (LxØ): 75 x 12 mm, with print**

**Page 1**



**Product description**

Order number	36.339
Product description	Sample tube, nominal volume: 2 ml, cap yellow, preparation: Fluoride + heparin, push cap, (LxØ): 75 x 12 mm, round base, with print, 100 piece(s)/bag

**Product characteristics**

Type of preparation	Fluoride + heparin
Preparation concentration	1.0 mg/ml blood + 16 IU/ml blood
Colour code	EU
Label/ Print	with print
Colour of print/label	yellow
Graduation	yes
Cap	Skirted push cap
Closure type	push cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Sample tube, Fluoride/heparin FH, 2 ml, cap yellow, (LxØ): 75 x 12 mm, with print**

**Page 2**

**Size**

Sample volume	2 ml
Diameter	12 mm
Length of product	75 mm
Length excluding cap	75 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	yellow
Tube material	Polypropylene (PP)

**Purity & certification**

Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	2000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	60000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0033 kg
Weight of case	7.31 kg
EAN of inner box	4038917081245
EAN case	4038917000383

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-06-30

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**