



### Product description

|                     |   |
|---------------------|---|
| Order number        | 41.1350.105   |
| Product description | Micro sample tube Citrate 9NC 0.106 mol/l 3.2%, venous blood collection, preparation: Citrate 9NC, nominal volume: 1.3 ml, (LxØ) with cap: 47 x 10.8 mm, screw cap, cap: blue, colour code ISO, with paper label, label/print: blue, skirted conical base, 100 piece(s)/bag |

### Product characteristics

|                           |                         |
|---------------------------|-------------------------|
| Type of collection        | venous                  |
| Type of preparation       | Citrate 9NC             |
| Preparation concentration | 0.106 mol/l             |
| Colour code               | ISO                     |
| Label/ Print              | with paper label        |
| Colour of print/label     | blue                    |
| Graduation                | yes                     |
| Closure type              | screw cap               |
| Application               | venous blood collection |
| Base shape                | skirted conical base    |
| Volume of work            | 1,3 ml                  |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**Micro sample tube Citrate 9NC 0.106 mol/l 3.2%,  
1.3 ml, screw cap, ISO**

**Page 2**

Recommended centrifugation 1500 x g  
Centrifugation time 10 Min

**Size**

Sample volume 1.3 ml  
Diameter 10.8 mm  
Length of product 47 mm  
Length including cap 47 mm  
Length excluding cap 45 mm

**Material & colours**

Colour of product transparent  
Closure material Polypropylene (PP)  
Colour of cap blue  
Tube material Polypropylene (PP)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

**Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 96000  
Depth of case 348 mm  
Width of case 278 mm  
Height of case 140 mm  
Case volume 0.0135 cbm  
Weight of product 0.0017 kg  
Weight of case 2.02 kg  
EAN of inner box 4038917127905  
EAN case 4038917021470

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**