



### Product description

Order number	41.1504.105
Product description	Micro sample tube EDTA K3E, venous blood collection, preparation: K3 EDTA, nominal volume: 1.3 ml, (LxØ) with cap: 44 x 10.8 mm, push cap, cap: violet, colour code ISO, with plastic label, label/print: violet, skirted conical base, 100 piece(s)/bag

### Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	violet
Graduation	yes
Cap	flat cap
Closure type	push cap
Application	venous blood collection
Base shape	skirted conical base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Volume of work 1,3 ml

### **Size**

Sample volume 1.3 ml  
Diameter 10.8 mm  
Length of product 44 mm  
Length including cap 44 mm  
Length excluding cap 40 mm

### **Material & colours**

Colour of product transparent  
Closure material Low Density Polyethylene (LD-PE)  
Colour of cap violet  
Tube material Polypropylene (PP)

### **Purity & certification**

Product category In vitrodiagnostic device, CE  
CE certified CE  
Batched yes

### **Packaging**

Minimum order qty. 1000  
Type of smallest subpackaging bag  
Piece(s) / inner box 100  
Piece (s) / outer case 1000  
Piece(s) / pallet 96000  
Depth of case 348 mm  
Width of case 278 mm  
Height of case 140 mm  
Case volume 0.0135 cbm  
Weight of product 0.0015 kg  
Weight of case 1.85 kg  
EAN of inner box 4038917129572  
EAN case 4038917028660

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

