



**Product Specification**  
**Sample tube, Serum CAT, 10 ml, cap white, (LxØ):**  
**101 x 16.5 mm**

**Page 1**



**Product description**

Order number	46.361
Product description	Sample tube, nominal volume: 10 ml, cap white, preparation: Clotting activator, screw cap, (LxØ): 101 x 16.5 mm, round base, 100 piece(s)/bag

**Product characteristics**

Type of preparation	Clotting activator
Colour code	EU
Graduation	no
Closure type	screw cap
Base shape	round base
Recommended centrifugation	2000 x g
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**Sample tube, Serum CAT, 10 ml, cap white, (LxØ):**  
**101 x 16.5 mm**

**Size**

Sample volume	10 ml
Diameter	16.5 mm
Length of product	101 mm
Length excluding cap	101 mm

**Material & colours**

Colour of product	transparent
Closure material	Low Density Polyethylene (LD-PE)
Colour of cap	white
Tube material	Polypropylene (PP)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	30000
Depth of case	398 mm
Width of case	398 mm
Height of case	340 mm
Case volume	0.0539 cbm
Weight of product	0.0069 kg
Weight of case	7.58 kg
EAN of inner box	4038917081313
EAN case	4038917000451

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

