



Product Specification
Sample tube, 10 ml, cap natural, (LxØ): 98.25 x 16.15 mm, with paper label



Product description

Order number	46.9924.106
Product description	Sample tube, nominal volume: 10 ml, PP, cap natural, screw cap, (LxØ): 98.25 x 16.15 mm, flat base, with paper label, 100 piece(s)/bag

Product characteristics

Label/ Print	with paper label
Colour of print/label	black
Graduation	no
Cap	standard screw cap
Closure type	screw cap
Application	collection and storage of diagnostic specimens
Base shape	flat base
Volume of work	10 ml
Centrifugation max (RCF)	5000 x g
Centrifugation temperature	21
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
Sample tube, 10 ml, cap natural, (LxØ): 98.25 x 16.15 mm, with paper label

Page 2

Size

Sample volume	10 ml
Diameter	16.15 mm
Length of product	98.25 mm

Material & colours

Product material	Polypropylene (PP)
Colour of product	translucent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	natural
Tube material	Polypropylene (PP)

Purity & certification

Satisfies the requirement	sealed at 18"hg/10min
Latex-free	yes
Pyrogen-free	no
Safety data sheet	yes

Packaging

Minimum order qty.	1000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	1000
Piece(s) / pallet	48000
Depth of case	604 mm
Width of case	264 mm
Height of case	373 mm
Case volume	0.06 cbm
Weight of case	7.6 kg
EAN of inner box	4038917154895
EAN case	4038917036467

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT