



**Product Specification**  
**Sample tube, EDTA K3E, 3 ml, cap violet, (LxØ):  
82 x 11.5 mm, with paper label**



**Product description**

Order number	47.556.100
Product description	Sample tube, nominal volume: 3 ml, cap violet, preparation: K3 EDTA, screw cap, (LxØ): 82 x 11.5 mm, skirted conical base, with paper label, 100 piece(s)/bag

**Product characteristics**

Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	violet
Graduation	fill mark
Closure type	screw cap
Base shape	skirted conical base

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



**Product Specification**  
**Sample tube, EDTA K3E, 3 ml, cap violet, (LxØ):**  
**82 x 11.5 mm, with paper label**

**Page 2**

**Size**

Sample volume	3 ml
Diameter	11.5 mm
Length of product	82 mm
Length excluding cap	66 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	violet
Tube material	Polypropylene (PP)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Batched	yes

**Packaging**

Minimum order qty.	50000
Type of smallest subpackaging	bag
Piece(s) / inner box	100
Piece (s) / outer case	2000
Piece(s) / pallet	72000
Depth of case	398 mm
Width of case	398 mm
Height of case	300 mm
Case volume	0.0475 cbm
Weight of product	0.0028 kg
Weight of case	6.37 kg
EAN of inner box	4038917113304
EAN case	4038917013987

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**