



**Product Specification**  
**S-Monovette® Serum CAT, 7.5 ml, cap white,**  
**(LxØ): 92 x 15 mm, with paper label**

**Page 1**



### Product description

Order number	01.1601.014
Product description	S-Monovette® Serum CAT, preparation: Clotting activator, 7.5 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 92 x 15 mm, with paper label "Blutbank", 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with paper label
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Serum CAT, 7.5 ml, cap white,**  
**(LxØ): 92 x 15 mm, with paper label**

**Page 2**

**Size**

Sample volume	7.5 ml
Outer diameter of cap	18 mm
Diameter	15 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	white
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitrodiagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	32000
Depth of box	189 mm
Width of box	134 mm
Height of box	82 mm
Depth of case	440 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Serum CAT, 7.5 ml, cap white,**  
**(LxØ): 92 x 15 mm, with paper label**

**Page 3**

Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0084 kg
Weight of case	5.04 kg
EAN of inner box	4038917131667
EAN case	4038917023849

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**