



**Product Specification**  
**S-Monovette® Serum CAT, 7.5 ml, cap red, (LxØ):**  
**92 x 15 mm, with paper label**

**Page 1**



### Product description

|                     |   |
|---------------------|---|
| Order number        | 01.1601.100   |
| Product description | S-Monovette® Serum CAT, preparation: Clotting activator, 7.5 ml, membrane screw cap, cap red, colour code ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: red, 50 piece(s)/case, sterile |

### Product characteristics

|                            |                    |
|----------------------------|--------------------|
| Type of collection         | venous             |
| Type of preparation        | Clotting activator |
| Colour code                | ISO                |
| Label/ Print               | with paper label   |
| Colour of print/label      | red                |
| Graduation                 | fill mark          |
| Cap                        | membrane screw cap |
| Closure type               | screw cap          |
| Recommended centrifugation | 2000 x g           |
| Centrifugation temperature | 20                 |
| Centrifugation time        | 10 Min             |
| Liquid-tight               | yes                |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Serum CAT, 7.5 ml, cap red, (LxØ):**  
**92 x 15 mm, with paper label**

**Page 2**

Connection membrane screw cap

**Size**

Sample volume 7.5 ml  
Outer diameter of cap 18 mm  
Diameter 15 mm  
Length including cap 108 mm  
Length excluding cap 92 mm  
Length including cap and plunger 129 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap red  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 32000  
Depth of box 189 mm  
Width of box 134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Serum CAT, 7.5 ml, cap red, (LxØ):**  
**92 x 15 mm, with paper label**

**Page 3**

|                   |               |
|-------------------|---------------|
| Height of box     | 82 mm         |
| Depth of case     | 440 mm        |
| Width of case     | 400 mm        |
| Height of case    | 144 mm        |
| Case volume       | 0.0253 cbm    |
| Weight of product | 0.0084 kg     |
| Weight of case    | 5.05 kg       |
| EAN of inner box  | 4038917131773 |
| EAN case          | 4038917023955 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**