



### Product description

Order number	01.1602.001
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 7.5 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 92 x 15 mm, with paper label, label/print: brown, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with paper label
Colour of print/label	brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Serum Gel CAT, 7.5 ml, cap brown,**  
**(LxØ): 92 x 15 mm, with paper label**

Connection membrane screw cap

**Size**

Sample volume 7.5 ml  
Outer diameter of cap 18 mm  
Diameter 15 mm  
Length including cap 108 mm  
Length excluding cap 92 mm  
Length including cap and plunger 129 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap brown  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 24000  
Depth of box 189 mm  
Width of box 134 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Serum Gel CAT, 7.5 ml, cap brown,  
(LxØ): 92 x 15 mm, with paper label**

**Page 3**

Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0088 kg
Weight of case	5.24 kg
EAN of inner box	4038917131797
EAN case	4038917023979

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**