



## Product Specification

**S-Monovette® Ammonium heparin AH, 7.5 ml, cap blue, (LxØ): 92 x 15 mm, with plastic label**

**Page 1**



### Product description

Order number	01.1603
Product description	S-Monovette® Ammonium heparin AH, preparation: Ammonium heparin, 7.5 ml, membrane screw cap, cap blue, colour code EU, (LxØ) without cap: 92 x 15 mm, with plastic label, label/print: transparent/blue, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Ammonium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Ammonium heparin AH, 7.5 ml, cap blue, (LxØ): 92 x 15 mm, with plastic label**

**Page 2**

Liquid-tight Connection yes  
membrane screw cap

### Size

Sample volume 7.5 ml  
Outer diameter of cap 18 mm  
Diameter 15 mm  
Length including cap 108 mm  
Length excluding cap 92 mm  
Length including cap and plunger 129 mm

### Material & colours

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap blue  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### Packaging

Minimum order qty. 500  
Type of smallest subpackaging case  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 32000  
Depth of box 189 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Ammonium heparin AH, 7.5 ml, cap blue, (LxØ): 92 x 15 mm, with plastic label**

**Page 3**

Width of box	134 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0084 kg
Weight of case	5.03 kg
EAN of inner box	4038917087216
EAN case	4038917005685

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**