



Product description

Order number	01.1604
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 7.5 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 92 x 15 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 7.5 ml, cap orange, (LxØ): 92 x 15 mm, with plastic label

Width of box	134 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	144 mm
Case volume	0.0253 cbm
Weight of product	0.0084 kg
Weight of case	5.03 kg
EAN of inner box	4038917087223
EAN case	4038917005692

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

