



## Product Specification

**S-Monovette® RNA Exact, max. 2.4 ml, cap green,  
(LxØ): 100 x 15 mm, with paper label**

**Page 1**



### Product description

Order number	01.2048.001
Product description	S-Monovette® RNA Exact, preparation: RNA stabiliser, max. 2.4 ml, membrane screw cap, cap green, (LxØ) without cap: 100 x 15 mm, with paper label, label/print: green, 20 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	RNA stabiliser
Label/ Print	with paper label
Colour of print/label	green
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® RNA Exact, max. 2.4 ml, cap green,  
(LxØ): 100 x 15 mm, with paper label**

**Page 2**

### Size

Max. volume	2.4 ml
Outer diameter of cap	18 mm
Diameter	15 mm
Length including cap	116 mm
Length excluding cap	100 mm
Length including cap and plunger	136 mm

### Material & colours

Product material	Polypropylene (PP)
Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	20
Type of smallest subpackaging	case
Piece(s) / inner box	20
Piece (s) / outer case	80
Piece(s) / pallet	8640
Depth of box	184 mm
Width of box	141 mm
Height of box	83 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® RNA Exact, max. 2.4 ml, cap green,  
(LxØ): 100 x 15 mm, with paper label**

**Page 3**

Depth of case	402 mm
Width of case	297 mm
Height of case	100 mm
Case volume	0.0119 cbm
Weight of product	0.0149 kg
Weight of case	1.86 kg
EAN of inner box	4038917550475
EAN case	4038917550468

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-16

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**