



Product Specification
**S-Monovette® EDTA K3E, 9 ml, cap violet, (LxØ):
92 x 16 mm, with paper label**



Product description

Order number	02.1066.100
Product description	S-Monovette® EDTA K3E, K3E, preparation: K3 EDTA, 9 ml, membrane screw cap, cap violet, colour code ISO, (LxØ) without cap: 92 x 16 mm, with paper label, label/print: violet, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	violet
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Centrifugation temperature	20
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® EDTA K3E, 9 ml, cap violet, (LxØ):
92 x 16 mm, with paper label

Page 2

Size

Sample volume	9 ml
Outer diameter of cap	19 mm
Diameter	16 mm
Length including cap	108 mm
Length excluding cap	92 mm
Length including cap and plunger	129 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	violet
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	24000
Depth of box	204 mm
Width of box	134 mm
Height of box	88 mm
Depth of case	648 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® EDTA K3E, 9 ml, cap violet, (LxØ):
92 x 16 mm, with paper label

Page 3

Width of case	303 mm
Height of case	150 mm
Case volume	0.0295 cbm
Weight of product	0.0079 kg
Weight of case	4.88 kg
EAN of inner box	4038917350334
EAN case	4038917350327

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT