



**Product Specification**  
**S-Monovette® EDTA K3E, 4 ml, cap red, (LxØ): 75 x 15 mm, with paper label**



**Product description**

Order number	03.1068.001
Product description	S-Monovette® EDTA K3E, preparation: K3 EDTA, 4 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: pink, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	K3 EDTA
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	pink
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® EDTA K3E, 4 ml, cap red, (LxØ): 75  
x 15 mm, with paper label**

**Page 2**

**Size**

Sample volume	4 ml
Outer diameter of cap	18 mm
Diameter	15 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	red
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

**Packaging**

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	34000
Depth of box	189 mm
Width of box	117 mm
Height of box	82 mm
Depth of case	440 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® EDTA K3E, 4 ml, cap red, (LxØ): 75  
x 15 mm, with paper label**

Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0062 kg
Weight of case	3.93 kg
EAN of inner box	4038917125253
EAN case	4038917019866

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

