



### Product description

Order number	03.1628
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 5.5 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 15 mm, with plastic label, label/print: transparent/orange, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/orange
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Lithium heparin LH, 5.5 ml, cap orange, (LxØ): 75 x 15 mm, with plastic label**

Liquid-tight yes  
Connection membrane screw cap

**Size**

Sample volume 5.5 ml  
Outer diameter of cap 18 mm  
Diameter 15 mm  
Length including cap 91 mm  
Length excluding cap 75 mm  
Length including cap and plunger 112 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap orange  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 34000  
Depth of box 189 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Lithium heparin LH, 5.5 ml, cap orange, (LxØ): 75 x 15 mm, with plastic label**

Width of box	117 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0069 kg
Weight of case	4.27 kg
EAN of inner box	4038917087285
EAN case	4038917005746

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

