




Product Specification

S-Monovette® Lithium heparin gel LH, 4.7 ml, cap green, (LxØ): 75 x 15 mm, with paper label

Page 1



Product description

Order number	03.1631.100
Product description	S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 4.7 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 75 x 15 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel
Preparation concentration	25 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Lithium heparin gel LH, 4.7 ml, cap green, (LxØ): 75 x 15 mm, with paper label

Page 3

Width of box	117 mm
Height of box	82 mm
Depth of case	440 mm
Width of case	400 mm
Height of case	129 mm
Case volume	0.0227 cbm
Weight of product	0.0075 kg
Weight of case	4.54 kg
EAN of inner box	4038917132213
EAN case	4038917024327

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT