



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 2.9 ml, cap green, (LxØ): 65 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix

Page 1



Product description

Order number	04.1902.088
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 2.9 ml, membrane screw cap, cap green, colour code EU, (LxØ) without cap: 65 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 015, label/print: white/green, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with plastic label pre-barcoded
Colour of print/label	white/green
Graduation	fill mark
Cap	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 2.9 ml, cap green, (LxØ): 65 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix

Page 2

Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	2.9 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	81 mm
Length excluding cap	65 mm
Length including cap and plunger	102 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	green
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 2.9 ml, cap green, (LxØ): 65 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix

Page 3

Packaging

Minimum order qty.	20000
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	42000
Depth of box	171 mm
Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.0042 kg
Weight of case	2.82 kg
EAN of inner box	4038917949446
EAN case	4038917949453

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT