



Product Specification
S-Monovette® Sodium heparin NH, 2.6 ml, cap green, (LxØ): 65 x 13 mm, with paper label

Page 1



Product description

Order number	04.1913.100
Product description	S-Monovette® Sodium heparin NH, preparation: Sodium heparin, 2.6 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 65 x 13 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Sodium heparin
Preparation concentration	20 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Sodium heparin NH, 2.6 ml, cap green, (LxØ): 65 x 13 mm, with paper label

Page 3

Width of box	105 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	125 mm
Case volume	0.0174 cbm
Weight of product	0.0043 kg
Weight of case	2.86 kg
EAN of inner box	4038917134330
EAN case	4038917072533

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT