




## Product Specification

**S-Monovette® Fluoride/EDTA FE, 2.7 ml, cap yellow, (LxØ): 75 x 13 mm, with plastic label**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 04.1918  |
| Product description | S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 2.7 ml, membrane screw cap, cap yellow, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/yellow, 50 piece(s)/case, sterile |

### Product characteristics

|                            |                                   |
|----------------------------|-----------------------------------|
| Type of collection         | venous                            |
| Type of preparation        | Fluoride + EDTA                   |
| Preparation concentration  | 1.0 mg/ml blood + 1.2 mg/ml blood |
| Colour code                | EU                                |
| Label/ Print               | with plastic label                |
| Colour of print/label      | transparent/yellow                |
| Graduation                 | fill mark                         |
| Cap                        | membrane screw cap                |
| Closure type               | screw cap                         |
| Recommended centrifugation | 2000 x g                          |
| Centrifugation temperature | 20                                |
| Centrifugation time        | 10 Min                            |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Fluoride/EDTA FE, 2.7 ml, cap yellow, (LxØ): 75 x 13 mm, with plastic label**

**Page 2**

Liquid-tight Connection      yes  
membrane screw cap

### Size

Sample volume                      2.7 ml  
Outer diameter of cap            16 mm  
Diameter                              13 mm  
Length including cap            91 mm  
Length excluding cap           75 mm  
Length including cap and plunger      112 mm

### Material & colours

Colour of product                      transparent  
Closure material                    High Density Polyethylene (HD-PE)  
Colour of cap                              yellow  
Membrane material                  rubber  
Tube material                              Polypropylene (PP)  
Piston material                        High Density Polyethylene (HD-PE)  
Piston plunger material            Polystyrene (PS)

### Purity & certification

Satisfies the requirement        IATA, ADR  
Product category                    In vitro diagnostic device, CE  
CE certified                              CE  
Purity standard                        sterile  
Sterilisation                            Electron irradiation  
Batched                                    yes

### Packaging

Minimum order qty.                   500  
Type of smallest subpackaging      case  
Piece(s) / inner box                50  
Piece (s) / outer case                500  
Piece(s) / pallet                       39000  
Depth of box                            171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Fluoride/EDTA FE, 2.7 ml, cap yellow, (LxØ): 75 x 13 mm, with plastic label**

**Page 3**

|                   |               |
|-------------------|---------------|
| Width of box      | 116 mm        |
| Height of box     | 73 mm         |
| Depth of case     | 388 mm        |
| Width of case     | 358 mm        |
| Height of case    | 135 mm        |
| Case volume       | 0.0188 cbm    |
| Weight of product | 0.0042 kg     |
| Weight of case    | 2.86 kg       |
| EAN of inner box  | 4038917226233 |
| EAN case          | 4038917226226 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**