




Product Specification
S-Monovette® Fluoride/EDTA FE, 2.7 ml, cap yellow, (LxØ): 75 x 13 mm, with paper label

Page 1



Product description

Order number	04.1918.001
Product description	S-Monovette® Fluoride/EDTA FE, preparation: Fluoride + EDTA, 2.7 ml, membrane screw cap, cap yellow, colour code EU, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: yellow, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Fluoride + EDTA
Preparation concentration	1.0 mg/ml blood + 1.2 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	yellow
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Fluoride/EDTA FE, 2.7 ml, cap yellow, (LxØ): 75 x 13 mm, with paper label

Page 3

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0043 kg
Weight of case	2.88 kg
EAN of inner box	4038917225144
EAN case	4038917225137

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT