



**Product Specification**  
**S-Monovette® EDTA K3E, 4 ml, cap violet, (LxØ):  
75 x 13 mm, with paper label**

**Page 1**



### Product description

|                     |  |
|---------------------|--|
| Order number        | 04.1921.100  |
| Product description | S-Monovette® EDTA K3E, K3E, preparation: K3 EDTA, 4 ml, membrane screw cap, cap violet, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: violet, 50 piece(s)/case, sterile |

### Product characteristics

|                           |                    |
|---------------------------|--------------------|
| Type of collection        | venous             |
| Type of preparation       | K3 EDTA            |
| Preparation concentration | 1.6 mg/ml blood    |
| Colour code               | ISO                |
| Label/ Print              | with paper label   |
| Colour of print/label     | violet             |
| Graduation                | fill mark          |
| Cap                       | membrane screw cap |
| Closure type              | screw cap          |
| Liquid-tight              | yes                |
| Connection                | membrane screw cap |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® EDTA K3E, 4 ml, cap violet, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 2**

**Size**

|                                  |        |
|----------------------------------|--------|
| Sample volume                    | 4 ml   |
| Outer diameter of cap            | 16 mm  |
| Diameter                         | 13 mm  |
| Length including cap             | 91 mm  |
| Length excluding cap             | 75 mm  |
| Length including cap and plunger | 112 mm |

**Material & colours**

|                         |                                   |
|-------------------------|-----------------------------------|
| Colour of product       | transparent                       |
| Closure material        | High Density Polyethylene (HD-PE) |
| Colour of cap           | violet                            |
| Membrane material       | rubber                            |
| Tube material           | Polypropylene (PP)                |
| Piston material         | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS)                  |

**Purity & certification**

|                           |                                |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR                      |
| Product category          | In vitro diagnostic device, CE |
| CE certified              | CE                             |
| Purity standard           | sterile                        |
| Sterilisation             | Electron irradiation           |
| Batched                   | yes                            |

**Packaging**

|                               |        |
|-------------------------------|--------|
| Minimum order qty.            | 20000  |
| Type of smallest subpackaging | case   |
| Piece(s) / inner box          | 50     |
| Piece (s) / outer case        | 500    |
| Piece(s) / pallet             | 39000  |
| Depth of box                  | 171 mm |
| Width of box                  | 116 mm |
| Height of box                 | 73 mm  |
| Depth of case                 | 388 mm |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® EDTA K3E, 4 ml, cap violet, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 3**

|                   |               |
|-------------------|---------------|
| Width of case     | 358 mm        |
| Height of case    | 135 mm        |
| Case volume       | 0.0188 cbm    |
| Weight of product | 0.0046 kg     |
| Weight of case    | 3.05 kg       |
| EAN of inner box  | 4038917162708 |
| EAN case          | 4038917162692 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**