



**Product Specification**  
**S-Monovette® Serum Gel CAT, 4 ml, cap brown,**  
**(LxØ): 75 x 13 mm, with plastic label**



**Product description**

Order number	04.1925
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 4 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label
Colour of print/label	transparent/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Serum Gel CAT, 4 ml, cap brown,**  
**(LxØ): 75 x 13 mm, with plastic label**

**Page 2**

Liquid-tight yes  
Connection membrane screw cap

**Size**

Sample volume 4 ml  
Outer diameter of cap 16 mm  
Diameter 13 mm  
Length including cap 91 mm  
Length excluding cap 75 mm  
Length including cap and plunger 112 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap brown  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 39000  
Depth of box 171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Serum Gel CAT, 4 ml, cap brown,**  
**(LxØ): 75 x 13 mm, with plastic label**

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0054 kg
Weight of case	3.47 kg
EAN of inner box	4038917248310
EAN case	4038917248303

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

