




**Product Specification**  
**S-Monovette® neutral Z, 4.9 ml, cap white, (LxØ):  
90 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	04.1926.001
Product description	S-Monovette® neutral Z, without preparation, 4.9 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: black, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	without preparation
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	black
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® neutral Z, 4.9 ml, cap white, (LxØ):**  
**90 x 13 mm, with paper label**

**Page 2**

Connection membrane screw cap

**Size**

Sample volume 4.9 ml  
Outer diameter of cap 16 mm  
Diameter 13 mm  
Length including cap 106 mm  
Length excluding cap 90 mm  
Length including cap and plunger 127 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap white  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 36000  
Depth of box 171 mm  
Width of box 130 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® neutral Z, 4.9 ml, cap white, (LxØ):**  
**90 x 13 mm, with paper label**

**Page 3**

Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0048 kg
Weight of case	3.18 kg
EAN of inner box	4038917134422
EAN case	4038917025683

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**