



## Product Specification

**S-Monovette® Lithium heparin gel LH, 4 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label**

**Page 1**



### Product description

Order number	04.1927
Product description	S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 4 ml, membrane screw cap, cap orange, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label, label/print: transparent/brown, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel
Preparation concentration	25 IU/ml blood
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/brown
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	15 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel LH, 4 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label**

Page 2

Liquid-tight yes  
Connection membrane screw cap

### Size

Sample volume 4 ml  
Outer diameter of cap 16 mm  
Diameter 13 mm  
Length including cap 91 mm  
Length excluding cap 75 mm  
Length including cap and plunger 112 mm

### Material & colours

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap orange  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material High Density Polyethylene (HD-PE)  
Piston plunger material Polystyrene (PS)

### Purity & certification

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

### Packaging

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 39000  
Depth of box 171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Lithium heparin gel LH, 4 ml, cap orange, (LxØ): 75 x 13 mm, with plastic label**

**Page 3**

Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0054 kg
Weight of case	3.45 kg
EAN of inner box	4038917394987
EAN case	4038917394970

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**