



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 5.4 ml, cap blue, (LxØ): 90 x 13 mm, with paper label**

**Page 1**

### Product description

Order number	04.1930.100
Product description	S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 5.4 ml, membrane screw cap, cap blue, colour code ISO, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: blue, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Citrate 9NC
Preparation concentration	0.106 mol/l
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	1800 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

### Size

Sample volume	5.4 ml
Outer diameter of cap	13 mm
Diameter	13 mm
Length including cap	106 mm
Length excluding cap	90 mm
Length including cap and plunger	127 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 5.4 ml, cap blue, (LxØ): 90 x 13 mm, with paper label**

Page 2

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	blue
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	500
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	36000
Depth of box	171 mm
Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0054 kg
Weight of case	3.46 kg
EAN of inner box	4038917964630
EAN case	4038917964647

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-11-13

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**