



Product Specification
S-Monovette® EDTA Gel K2E, 4.9 ml, cap red,
(LxØ): 90 x 13 mm, with paper label

Page 1



Product description

Order number	04.1932.001
Product description	S-Monovette® EDTA Gel K2E, preparation: K2 EDTA Gel, 4.9 ml, membrane screw cap, cap red, colour code EU, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: red, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	K2 EDTA Gel
Preparation concentration	1.6 mg/ml blood
Colour code	EU
Label/ Print	with paper label
Colour of print/label	red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® EDTA Gel K2E, 4.9 ml, cap red,
(LxØ): 90 x 13 mm, with paper label

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 4.9 ml
Outer diameter of cap 16 mm
Diameter 13 mm
Length including cap 106 mm
Length excluding cap 90 mm
Length including cap and plunger 127 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap red
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 36000
Depth of box 171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® EDTA Gel K2E, 4.9 ml, cap red,
(LxØ): 90 x 13 mm, with paper label

Page 3

Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0056 kg
Weight of case	3.58 kg
EAN of inner box	4038917134521
EAN case	4038917028868

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT