



Product Specification

**S-Monovette® GlucoEXACT FC, 3.1 ml, cap pink,
(LxØ): 75 x 13 mm, with paper label**

Page 1

Product description

Order number	04.1945.201
Product description	S-Monovette® GlucoEXACT FC, preparation: Citrate/fluoride/EDTA, 3.1 ml, membrane screw cap, cap pink, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: pink, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Citrate/fluoride/EDTA
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	pink
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

Size

Sample volume	3.1 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	pink
Membrane material	rubber
Tube material	Polypropylene (PP)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP
support and is valid without signature.



SARSTEDT



Product Specification

**S-Monovette® GlucoEXACT FC, 3.1 ml, cap pink,
(LxØ): 75 x 13 mm, with paper label**

Page 2

Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 20000
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 39000
Depth of box 171 mm
Width of box 116 mm
Height of box 73 mm
Depth of case 388 mm
Width of case 358 mm
Height of case 135 mm
Case volume 0.0188 cbm
Weight of product 0.0047 kg
Weight of case 3.08 kg
EAN of inner box 4038917554091
EAN case 4038917554084

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-17

This document was prepared by EDP support and is valid without signature.



SARSTEDT