



Product Specification
**S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
75 x 13 mm, with paper label**

Page 1

Product description

Order number 04.1947.791
Product description S-Monovette® neutral Z, without preparation, 2.7 ml, membrane screw cap, cap white, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: black, 50 piece(s)/case, sterile

Product characteristics

Type of collection venous
Type of preparation without preparation
Colour code ISO
Label/ Print with paper label
Colour of print/label black
Graduation fill mark
Cap membrane screw cap
Closure type screw cap
Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 2.7 ml
Outer diameter of cap 16 mm
Diameter 13 mm
Length including cap 91 mm
Length excluding cap 75 mm
Length including cap and plunger 112 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap white
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® neutral Z, 2.7 ml, cap white, (LxØ):
75 x 13 mm, with paper label

Page 2

Purity & certification

Satisfies the requirement	IATA, ADR
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

Packaging

Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Depth of box	170 mm
Width of box	115 mm
Height of box	72 mm
Depth of case	380 mm
Width of case	350 mm
Height of case	120 mm
Case volume	0.016 cbm
EAN of inner box	4038917959841
EAN case	4038917959858

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-11-13

This document was prepared by EDP support and is valid without signature.



SARSTEDT