



**Product Specification**  
**S-Monovette® DNA Exact, 2.7 ml, cap blue, (LxØ):**  
**75 x 13 mm, with paper label**



**Product description**

Order number	04.1948.001
Product description	S-Monovette® DNA Exact, Suitable for stabilising genomic DNA, preparation: EDTA based gDNA stabilizer, 2.7 ml, membrane screw cap, cap blue, colour code EU/ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: blue, venous blood collection, With validated time and temperature parameters, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	EDTA based gDNA stabilizer
Colour code	EU/ISO
Label/ Print	with paper label
Colour of print/label	blue
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Design	With validated time and temperature parameters
Application	venous blood collection

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.



**Product Specification**  
**S-Monovette® DNA Exact, 2.7 ml, cap blue, (LxØ):**  
**75 x 13 mm, with paper label**

**Page 2**

Example of use	Suitable for stabilising genomic DNA
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes
Connection	membrane screw cap

**Size**

Sample volume	2.7 ml
Outer diameter of cap	16 mm
Diameter	13 mm
Length including cap	91 mm
Length excluding cap	75 mm
Length including cap and plunger	112 mm

**Material & colours**

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	blue
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

**Purity & certification**

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® DNA Exact, 2.7 ml, cap blue, (LxØ):**  
**75 x 13 mm, with paper label**

**Packaging**

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	39000
Depth of box	171 mm
Width of box	116 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	135 mm
Case volume	0.0188 cbm
Weight of product	0.0042 kg
Weight of case	2.86 kg
EAN of inner box	4038917370745
EAN case	4038917370738

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

