



Product description

| | |
|---------------------|---|
| Order number | 04.1953.100 |
| Product description | S-Monovette® Lithium heparin gel+ LH, preparation: Lithium heparin gel+, 4 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 75 x 13 mm, with paper label, label/print: green, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|----------------------|
| Type of collection | venous |
| Type of preparation | Lithium heparin gel+ |
| Preparation concentration | 25 IU/ml blood |
| Colour code | ISO |
| Label/ Print | with paper label |
| Colour of print/label | green |
| Graduation | fill mark |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2500 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 7 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Liquid-tight Connection yes
membrane screw cap

Size

Sample volume 4 ml
Outer diameter of cap 16 mm
Diameter 13 mm
Length including cap 91 mm
Length excluding cap 75 mm
Length including cap and plunger 112 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap green
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest subpackaging case
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 39000
Depth of box 171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin gel+ LH, 4 ml, cap green, (LxØ): 75 x 13 mm, with paper label

| | |
|-------------------|---------------|
| Width of box | 116 mm |
| Height of box | 73 mm |
| Depth of case | 388 mm |
| Width of case | 358 mm |
| Height of case | 135 mm |
| Case volume | 0.0188 cbm |
| Weight of product | 0.0054 kg |
| Weight of case | 3.47 kg |
| EAN of inner box | 4038917523318 |
| EAN case | 4038917523301 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

