



## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 4.9 ml,  
cap green, (LxØ): 90 x 13 mm, with paper label**

**Page 1**



### Product description

Order number	04.1954.100
Product description	S-Monovette® Lithium heparin gel+ LH, preparation: Lithium heparin gel+, 4.9 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 90 x 13 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin gel+
Preparation concentration	25 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	7 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



# Product Specification

**S-Monovette® Lithium heparin gel+ LH, 4.9 ml,  
cap green, (LxØ): 90 x 13 mm, with paper label**

Liquid-tight Connection                                yes  
membrane screw cap

**Size**

Sample volume     4.9 ml  
Outer diameter of cap                                   16 mm  
Diameter    13 mm  
Length including cap                                    106 mm  
Length excluding cap                                    90 mm  
Length including cap and plunger                   127 mm

**Material & colours**

Colour of product                                        transparent  
Closure material                                        High Density Polyethylene (HD-PE)  
Colour of cap     green  
Membrane material                                      rubber  
Tube material    Polypropylene (PP)  
Piston material    High Density Polyethylene (HD-PE)  
Piston plunger material                                Polystyrene (PS)

**Purity & certification**

Satisfies the requirement                            IATA, ADR  
Product category                                        In vitrodiagnostic device, CE  
CE certified    CE  
Purity standard    sterile  
Sterilisation     Electron irradiation  
Batched    yes

**Packaging**

Minimum order qty.                                    500  
Type of smallest subpackaging                    case  
Piece(s) / inner box                                  50  
Piece (s) / outer case                                500  
Piece(s) / pallet                                        36000  
Depth of box    171 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





## Product Specification

**S-Monovette® Lithium heparin gel+ LH, 4.9 ml,  
cap green, (LxØ): 90 x 13 mm, with paper label**

**Page 3**

Width of box	130 mm
Height of box	73 mm
Depth of case	388 mm
Width of case	358 mm
Height of case	150 mm
Case volume	0.0208 cbm
Weight of product	0.0057 kg
Weight of case	3.63 kg
EAN of inner box	4038917522793
EAN case	4038917522786

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-03

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**