



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.8 ml, cap green, (LxØ): 75 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

Page 1



### Product description

|                     |  |
|---------------------|--|
| Order number        | 04.1955.088  |
| Product description | S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 1.8 ml, membrane screw cap, cap green, colour code EU, (LxØ) without cap: 75 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix 019, label/print: white/green, 50 piece(s)/case, sterile |

### Product characteristics

|                           |                                 |
|---------------------------|---------------------------------|
| Type of collection        | venous                          |
| Type of preparation       | Citrate 9NC                     |
| Preparation concentration | 0.106 mol/l                     |
| Colour code               | EU                              |
| Label/ Print              | with plastic label pre-barcoded |
| Colour of print/label     | white/green                     |
| Graduation                | fill mark                       |
| Cap                       | membrane screw cap              |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.8 ml, cap green, (LxØ): 75 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

Page 2

|                            |                    |
|----------------------------|--------------------|
| Closure type               | screw cap          |
| Recommended centrifugation | 3000 x g           |
| Centrifugation temperature | 20                 |
| Centrifugation time        | 10 Min             |
| Liquid-tight               | yes                |
| Connection                 | membrane screw cap |

### Size

|                                  |        |
|----------------------------------|--------|
| Sample volume                    | 1.8 ml |
| Outer diameter of cap            | 16 mm  |
| Diameter                         | 13 mm  |
| Length including cap             | 91 mm  |
| Length excluding cap             | 75 mm  |
| Length including cap and plunger | 112 mm |

### Material & colours

|                         |                                   |
|-------------------------|-----------------------------------|
| Colour of product       | transparent                       |
| Closure material        | High Density Polyethylene (HD-PE) |
| Colour of cap           | green                             |
| Membrane material       | rubber                            |
| Tube material           | Polypropylene (PP)                |
| Piston material         | High Density Polyethylene (HD-PE) |
| Piston plunger material | Polystyrene (PS)                  |

### Purity & certification

|                           |                                |
|---------------------------|--------------------------------|
| Satisfies the requirement | IATA, ADR                      |
| Product category          | In vitro diagnostic device, CE |
| CE certified              | CE                             |
| Purity standard           | sterile                        |
| Sterilisation             | Electron irradiation           |
| Batched                   | yes                            |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.8 ml, cap green, (LxØ): 75 x 13 mm, with plastic label pre-barcoded, Pre-barcode with 8-digit unique number range and 3-digit prefix**

**Page 3**

### Packaging

|                               |               |
|-------------------------------|---------------|
| Minimum order qty.            | 20000         |
| Type of smallest subpackaging | case          |
| Piece(s) / inner box          | 50            |
| Piece (s) / outer case        | 500           |
| Piece(s) / pallet             | 39000         |
| Depth of box                  | 171 mm        |
| Width of box                  | 116 mm        |
| Height of box                 | 73 mm         |
| Depth of case                 | 388 mm        |
| Width of case                 | 358 mm        |
| Height of case                | 135 mm        |
| Case volume                   | 0.0188 cbm    |
| Weight of product             | 0.0048 kg     |
| Weight of case                | 3.16 kg       |
| EAN of inner box              | 4038917949507 |
| EAN case                      | 4038917949514 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-18

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**