



## Product Specification

**S-Monovette® ThromboExact, 2.7 ml, cap dark red, (LxØ): 66 x 11 mm, with paper label**

**Page 1**



### Product description

Order number	05.1168.001
Product description	S-Monovette® ThromboExact, preparation: Mg2+ compound, 2.7 ml, membrane screw cap, cap dark red, (LxØ) without cap: 66 x 11 mm, with paper label, label/print: dark red, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Mg2+ compound
Preparation concentration	0.82 mg/ml blood
Label/ Print	with paper label
Colour of print/label	dark red
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Liquid-tight	yes
Connection	membrane screw cap

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® ThromboExact, 2.7 ml, cap dark red, (LxØ): 66 x 11 mm, with paper label**

**Page 2**

### Size

Sample volume	2.7 ml
Outer diameter of cap	13 mm
Diameter	11 mm
Length including cap	82 mm
Length excluding cap	66 mm
Length including cap and plunger	103 mm

### Material & colours

Colour of product	transparent
Closure material	High Density Polyethylene (HD-PE)
Colour of cap	dark red
Membrane material	rubber
Tube material	Polypropylene (PP)
Piston material	High Density Polyethylene (HD-PE)
Piston plunger material	Polystyrene (PS)

### Purity & certification

Satisfies the requirement	IATA, ADR
Product category	In vitro diagnostic device, CE
CE certified	CE
Purity standard	sterile
Sterilisation	Electron irradiation
Batched	yes

### Packaging

Minimum order qty.	50
Type of smallest subpackaging	case
Piece(s) / inner box	50
Piece (s) / outer case	500
Piece(s) / pallet	56000
Depth of box	145 mm
Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



## Product Specification

**S-Monovette® ThromboExact, 2.7 ml, cap dark red, (LxØ): 66 x 11 mm, with paper label**

**Page 3**

Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0039 kg
Weight of case	2.48 kg
EAN of inner box	4038917159906
EAN case	4038917159890

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-17

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**