



Product Specification
S-Monovette® Lithium heparin LH, 2.7 ml, cap green, (LxØ): 66 x 11 mm, with paper label



Product description

Order number	05.1553.100
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 2.7 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 66 x 11 mm, with paper label, label/print: green, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	16 IU/ml blood
Colour code	ISO
Label/ Print	with paper label
Colour of print/label	green
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Lithium heparin LH, 2.7 ml, cap green, (LxØ): 66 x 11 mm, with paper label

Page 3

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0042 kg
Weight of case	2.64 kg
EAN of inner box	4038917130431
EAN case	4038917031912

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-01

This document was prepared by EDP support and is valid without signature.



SARSTEDT