




Product Specification
S-Monovette® Serum CAT, 2.7 ml, cap white,
(LxØ): 66 x 11 mm, with plastic label

Page 1



Product description

Order number	05.1557
Product description	S-Monovette® Serum CAT, preparation: Clotting activator, 2.7 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 66 x 11 mm, with plastic label, label/print: transparent/white, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	transparent/white
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Serum CAT, 2.7 ml, cap white,
(LxØ): 66 x 11 mm, with plastic label

Page 2

Connection membrane screw cap

Size

Sample volume 2.7 ml
Outer diameter of cap 13 mm
Diameter 11 mm
Length including cap 82 mm
Length excluding cap 66 mm
Length including cap and plunger 103 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap white
Membrane material rubber
Tube material Polypropylene (PP)
Piston material High Density Polyethylene (HD-PE)
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 56000
Depth of box 145 mm
Width of box 108 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Product Specification
S-Monovette® Serum CAT, 2.7 ml, cap white,
(LxØ): 66 x 11 mm, with plastic label

Page 3

Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0041 kg
Weight of case	2.61 kg
EAN of inner box	4038917087025
EAN case	4038917005517

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT