



**Product Specification**  
**S-Monovette® Serum CAT, 1.2 ml, cap white,**  
**(LxØ): 66 x 8 mm, with plastic label**



### Product description

Order number	06.1663.001
Product description	S-Monovette® Serum CAT, preparation: Clotting activator, 1.2 ml, membrane screw cap, cap white, colour code EU, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/black, 50 piece(s)/case, sterile

### Product characteristics

Type of collection	venous
Type of preparation	Clotting activator
Colour code	EU
Label/ Print	with plastic label
Colour of print/label	white/black
Graduation	fill mark
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min
Liquid-tight	yes

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Serum CAT, 1.2 ml, cap white,**  
**(LxØ): 66 x 8 mm, with plastic label**

Connection membrane screw cap

**Size**

Sample volume 1.2 ml  
Outer diameter of cap 13 mm  
Diameter 8 mm  
Length including cap 82 mm  
Length excluding cap 66 mm  
Length including cap and plunger 103 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap white  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material rubber  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 56000  
Depth of box 145 mm  
Width of box 108 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Serum CAT, 1.2 ml, cap white,**  
**(LxØ): 66 x 8 mm, with plastic label**

Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0032 kg
Weight of case	2.16 kg
EAN of inner box	4038917132770
EAN case	4038917024419

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

