



**Product Specification**  
**S-Monovette® Lithium heparin LH, 1.2 ml, cap green, (LxØ): 66 x 8 mm, with plastic label**



**Product description**

Order number	06.1666.100
Product description	S-Monovette® Lithium heparin LH, preparation: Lithium heparin, 1.2 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/green, 50 piece(s)/case, sterile

**Product characteristics**

Type of collection	venous
Type of preparation	Lithium heparin
Preparation concentration	35 IU/ml blood
Colour code	ISO
Label/ Print	with plastic label
Colour of print/label	white/green
Graduation	yes
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2000 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





**Product Specification**  
**S-Monovette® Lithium heparin LH, 1.2 ml, cap green, (LxØ): 66 x 8 mm, with plastic label**

**Page 2**

Liquid-tight yes  
Connection membrane screw cap

**Size**

Sample volume 1.2 ml  
Outer diameter of cap 13 mm  
Diameter 8 mm  
Length including cap 82 mm  
Length excluding cap 66 mm  
Length including cap and plunger 103 mm

**Material & colours**

Colour of product transparent  
Closure material High Density Polyethylene (HD-PE)  
Colour of cap green  
Membrane material rubber  
Tube material Polypropylene (PP)  
Piston material rubber  
Piston plunger material Polystyrene (PS)

**Purity & certification**

Satisfies the requirement IATA, ADR  
Product category In vitrodiagnostic device, CE  
CE certified CE  
Purity standard sterile  
Sterilisation Electron irradiation  
Batched yes

**Packaging**

Minimum order qty. 500  
Type of smallest case  
subpackaging  
Piece(s) / inner box 50  
Piece (s) / outer case 500  
Piece(s) / pallet 56000  
Depth of box 145 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:  
2024-07-19

This document was prepared by EDP support and is valid without signature.



**SARSTEDT**



**Product Specification**  
**S-Monovette® Lithium heparin LH, 1.2 ml, cap green, (LxØ): 66 x 8 mm, with plastic label**

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0032 kg
Weight of case	2.14 kg
EAN of inner box	4038917132886
EAN case	4038917024525

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

