




Product Specification
S-Monovette® Serum Gel CAT, 1.1 ml, cap brown,
(LxØ): 66 x 8 mm, with plastic label

Page 1



Product description

Order number	06.1667.001
Product description	S-Monovette® Serum Gel CAT, preparation: Clotting activator/gel, 1.1 ml, membrane screw cap, cap brown, colour code EU/ISO, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/brown, 50 piece(s)/case, sterile

Product characteristics

Type of collection	venous
Type of preparation	Clotting activator/gel
Colour code	EU/ISO
Label/ Print	with plastic label
Colour of print/label	white/brown
Graduation	yes
Cap	membrane screw cap
Closure type	screw cap
Recommended centrifugation	2500 x g
Centrifugation temperature	20
Centrifugation time	10 Min

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Serum Gel CAT, 1.1 ml, cap brown,
(LxØ): 66 x 8 mm, with plastic label

Page 3

Width of box	108 mm
Height of box	63 mm
Depth of case	333 mm
Width of case	303 mm
Height of case	125 mm
Case volume	0.0126 cbm
Weight of product	0.0041 kg
Weight of case	2.66 kg
EAN of inner box	4038917132893
EAN case	4038917024532

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-18

This document was prepared by EDP support and is valid without signature.



SARSTEDT