



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.4 ml, cap green, (LxØ): 66 x 8 mm, with plastic label

Page 1



Product description

| | |
|---------------------|---|
| Order number | 06.1668.001 |
| Product description | S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, preparation: Citrate 9NC, 1.4 ml, membrane screw cap, cap green, colour code EU, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: green, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|--------------------|
| Type of collection | venous |
| Type of preparation | Citrate 9NC |
| Preparation concentration | 0.106 mol/l |
| Colour code | EU |
| Label/ Print | with plastic label |
| Colour of print/label | green |
| Graduation | yes |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 1800 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 10 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Citrate 9NC 0.106 mol/l 3.2%, 1.4 ml, cap green, (LxØ): 66 x 8 mm, with plastic label

Page 3

| | |
|-------------------|---------------|
| Width of box | 108 mm |
| Height of box | 63 mm |
| Depth of case | 333 mm |
| Width of case | 303 mm |
| Height of case | 125 mm |
| Case volume | 0.0126 cbm |
| Weight of product | 0.0033 kg |
| Weight of case | 2.2 kg |
| EAN of inner box | 4038917132909 |
| EAN case | 4038917024549 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-03

This document was prepared by EDP support and is valid without signature.



SARSTEDT