



Product Specification

S-Monovette® Lithium heparin gel LH, 1.1 ml, cap green, (LxØ): 66 x 8 mm, with plastic label

Page 1



Product description

| | |
|---------------------|--|
| Order number | 06.1669.100 |
| Product description | S-Monovette® Lithium heparin gel LH, preparation: Lithium heparin gel, 1.1 ml, membrane screw cap, cap green, colour code ISO, (LxØ) without cap: 66 x 8 mm, with plastic label, label/print: white/green, 50 piece(s)/case, sterile |

Product characteristics

| | |
|----------------------------|---------------------|
| Type of collection | venous |
| Type of preparation | Lithium heparin gel |
| Preparation concentration | 35 IU/ml blood |
| Colour code | ISO |
| Label/ Print | with plastic label |
| Colour of print/label | white/green |
| Graduation | yes |
| Cap | membrane screw cap |
| Closure type | screw cap |
| Recommended centrifugation | 2500 x g |
| Centrifugation temperature | 20 |
| Centrifugation time | 15 Min |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification

S-Monovette® Lithium heparin gel LH, 1.1 ml, cap green, (LxØ): 66 x 8 mm, with plastic label

Page 2

Liquid-tight yes
Connection membrane screw cap

Size

Sample volume 1.1 ml
Outer diameter of cap 13 mm
Diameter 8 mm
Length including cap 82 mm
Length excluding cap 66 mm
Length including cap and plunger 103 mm

Material & colours

Colour of product transparent
Closure material High Density Polyethylene (HD-PE)
Colour of cap green
Membrane material rubber
Tube material Polypropylene (PP)
Piston material rubber
Piston plunger material Polystyrene (PS)

Purity & certification

Satisfies the requirement IATA, ADR
Product category In vitrodiagnostic device, CE
CE certified CE
Purity standard sterile
Sterilisation Electron irradiation
Batched yes

Packaging

Minimum order qty. 500
Type of smallest case
subpackaging
Piece(s) / inner box 50
Piece (s) / outer case 500
Piece(s) / pallet 56000
Depth of box 145 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT



Product Specification
S-Monovette® Lithium heparin gel LH, 1.1 ml, cap green, (LxØ): 66 x 8 mm, with plastic label

Page 3

| | |
|-------------------|---------------|
| Width of box | 108 mm |
| Height of box | 63 mm |
| Depth of case | 333 mm |
| Width of case | 303 mm |
| Height of case | 125 mm |
| Case volume | 0.0126 cbm |
| Weight of product | 0.0033 kg |
| Weight of case | 2.19 kg |
| EAN of inner box | 4038917271837 |
| EAN case | 4038917271820 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

Date of issue:
2024-07-19

This document was prepared by EDP support and is valid without signature.



SARSTEDT