



Product description

| | |
|---------------------|--|
| Order number | 83.3950 |
| Product description | Cell scraper, 2-position blade, size: S, handle length: 24 cm, blade length: 1.35 cm, sterile, non-pyrogenic/endotoxin-free, non-cytotoxic, 1 piece(s)/blister |

Size

| | |
|-------------------|---------|
| Length of product | 24.7 cm |
| Length of blade | 1.35 cm |
| Length of stem | 24 cm |
| Width of blade | 1.35 cm |

Material & colours

| | |
|-------------------|--------------------------------|
| Product material | Polystyrene (PS) |
| Colour of product | red/opaque |
| Stem material | Polystyrene (PS) |
| Blade material | Thermoplastic elastomers (TPE) |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|--|
| Product category | no medical device IVD |
| Purity standard | sterile, non-pyrogenic/endotoxin-free, non-cytotoxic |
| Sterilisation | Electron irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 100 |
| Type of smallest subpackaging | blister |
| Piece(s) / inner box | 1 |
| Piece (s) / outer case | 100 |
| Piece(s) / pallet | 16500 |
| Depth of case | 276 mm |
| Width of case | 110 mm |
| Height of case | 333 mm |
| Case volume | 0.0101 cbm |
| Weight of product | 0.0081 kg |
| Weight of case | 1.3 kg |
| EAN of inner box | 4038917426633 |
| EAN case | 4038917426626 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

