



### **Product description**

Order number	94.6077.435
Product description	flexiPERM® , con B, re-usable silicon insert, for dividing the growth surface, 1 well, non-pyrogenic/endotoxin-free, non-cytotoxic, 5 piece(s)/bag

### **Product characteristics**

Application	for dividing the growth surface
Number of wells	1
Growth area	3.1 cm <sup>2</sup>
Volume of work	2 ml - 3 ml

### **Size**

Diameter	20 mm
Height of product	8 mm

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





### **Material & colours**

Product material silicone  
Colour of product transparent

### **Purity & certification**

Product category no medical device | IVD  
Purity standard non-pyrogenic/endotoxin-free, non-cytotoxic  
Non-cytotoxic yes  
Batched yes

### **Packaging**

Minimum order qty. 5  
Type of smallest subpackaging bag  
Piece(s) / inner box 5  
Piece (s) / outer case 100  
Piece(s) / pallet 7200  
Depth of case 388 mm  
Width of case 253 mm  
Height of case 225 mm  
Case volume 0.0221 cbm  
Weight of product 0.0039 kg  
Weight of case 1.42 kg  
EAN of inner box 4038917255295  
EAN case 4038917255288

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

