



Product description

| | |
|---------------------|---|
| Order number | 83.3912.002 |
| Product description | Cell culture flask, T-175, material: PS, surface: Standard, for Adherent cells, coding colour: red, Filter cap, TC Tested, 5 piece(s)/bag |

Product characteristics

| | |
|----------------|---------------------|
| Cap | Filter cap |
| Cell type | Adherent cells |
| Growth area | 175 cm ² |
| Volume of work | 50 ml - 125 ml |
| Vessel type | flask |
| Coding colour | red |

Material & colours

| | |
|------------------|-------------------|
| Product material | Polystyrene (PS) |
| Closure material | Polyethylene (PE) |
| Colour of cap | red |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.





Purity & certification

| | |
|------------------|-------------------------|
| Product category | no medical device IVD |
| Purity standard | TC Tested |
| Sterilisation | Gamma irradiation |
| Pyrogen-free | yes |
| Non-cytotoxic | yes |
| Batched | yes |

Packaging

| | |
|-------------------------------|---------------|
| Minimum order qty. | 40 |
| Type of smallest subpackaging | bag |
| Piece(s) / inner box | 5 |
| Piece (s) / outer case | 40 |
| Piece(s) / pallet | 960 |
| Depth of case | 522 mm |
| Width of case | 267 mm |
| Height of case | 450 mm |
| Case volume | 0.0627 cbm |
| Weight of product | 0.1057 kg |
| Weight of case | 5.62 kg |
| EAN of inner box | 4038917337748 |
| EAN case | 4038917337731 |

This is the current specification for this product. Sarstedt reserves the right to make changes, in full or in part, at any time without prior notification.

